

## EHBT-75 Software Update

Dear customers, we have released the new software version **V4.0.18** of EHBT-75, the main updated contents are as follows:

### **1. Testing:**

1.1 HCT coefficient added for dry type QC and liquid QC.

### **2. Report:**

2.1 The report export options have been updated to include three options: PDF, 'Spreadsheet' and 'Image Pack'.

2.2 Added language Russian.

### **3. Workbench:**

3.1 Receives and displays data returned from Open Dx review. And order deletion linked.

This version was released on 2026/5/7. Please connect to the network and update the version ASAP. If the network at your site is poor, please contact us and we will send you the software package and guide you on the upgrade.