

EHBT-25 Software Update

Dear customers, we have released the new software version **V4.0.19** of

EHBT-25, the main updated contents are as follows:

1. Report:

1.1 Optimization of clinic information settings: Removed the original "Remarks" field and added "Clinic Phone Number" and "Clinic Address" fields.

2. Print:

2.1 Thermal printing supports portrait orientation; landscape orientation is selected by default.

2.2 A4 supports minimalist printing mode; the standard A4 printing mode is selected by default.

2.3 Print settings: All options other than the clinic address and telephone number are ticked by default.

3. Others:

3.1 French is supported.

This version was released on 2026/5/26. Please connect to the network and update the version ASAP. If the network at your site is poor, please contact us and we will send you the software package and guide you on the upgrade.