

EHBT-50

Software Update Notification Version 4.0.21

V 4.0.21	2026-06-17	3 Updated	1
RELEASE	RELEASE DATE	MODULES	NEW FEATURES

Dear Customer,

To continuously improve our products and elevate your experience, we are pleased to release version **4.0.21** for the EHBT-50. The following are the main contents of this update.

Report

OPTIMIZED Optimization of RBC histogram

Quality Control (QC)

OPTIMIZED Optimization of CBC Dry QC


Settings

NEW Equipment maintenance

Added "Equipment maintenance" entry, password: 122145.

Added the "Biochemistry Sample Dispensing Calibration".

"Imaging system QC" entry has been moved to "Equipment maintenance".



This version has been officially released. You can connect your device to the network, wait for the software package to download, and after seeing the update prompt, proceed with the online upgrade.

If you are unable to upgrade online, you can download the software package from the following link:

https://drive.google.com/drive/folders/1cvPnEm9di1bULwt3NDzhtuL5y8-RIQkG?usp=drive_link

Pls contact the after-sales team to guide you to perform the upgrade via a USB drive.

Thank you for choosing Ozelle.

For any questions about this release, please contact the after-sale service team for support.