

EHVT-50 Software Update Points

Dear customers, we have released the new software version **V4.0.21** of EHVT-50, the main updated contents are as follows:

1. Further optimization of urine and fecal detection.
2. Report export: 1) Modify the naming rule for the exported data.
2) Separate export tables for fecal, urine, and blood test data to avoid confusion and improve data management.
3. Report Details: 1) Modification of the abbreviations for fecal test report (parasite eggs, intestinal protozoa, pathogenic microorganisms, and digestive products).
2) Add the fluorescence curve when exporting PDF reports with immune detection.
4. Consumables: 1) The name of the calibration has be changed to consumables.
2) Modification of relevant text on the consumables page.
5. Add "Report viewing address" function in Settings: Customers can view, download and print PDF reports directly on a computer.
6. The Portuguese translation of the "View" button in the report list has been modified to: "Ver".

This version was released on December 10, 2025. Please connect to the network and update the version ASAP. If the network at your site is poor, please contact us and we will send you the software package and guide you on the upgrade.